

	Week Beginning:					Week Beginning:					Week Beginning:					Week Beginning:					
	Mon	Tue	Wed	Thur	Fri	Mon	Tue	Wed	Thur	Fri	Mon	Tue	Wed	Thur	Fri	Mon	Tue	Wed	Thur	Fri	
Check sample path flow																					
Check Turbo Pump Frequency																					
Check Ion Gauge <2																					
Check Convector Gauge																					
Check Waste Levels																					
Check HPLC Solvents																					
Check Sampler wash solvent																					
Record Fridge Temperatures																					
Initial																					
Run Fresh QC																					
Run sensitivity Check																					
Check Pending List																					
Wipe Benches																					
Initial																					
Check Helium Pressure																					
Check Pump Oil Level																					
Release Pump Oil Pressure																					
Remove & Clean Sweep Cone																					
Clean Transfer Capillary																					
Record Column Backpressure																					
Check ESI Cal response																					
Clean Tube Lens & Skimmer																					
Backup HPLC & MS Methods																					
Clean Air Filter																					